



HAORUN MEDICAL PRODUCTS CO.,LTD
TEL 0086-574-87706193
WHATSAPP/WECHAT/TEL:008613566542343
EMAIL:SALES@HAORUNMED.COM
WWW.HAORUNMED.COM





Forna Medical Instrument

Infusion Set



Item code:JF001
 Specification :
 Air vented spike with cap
 48mm drip with $\text{Ø}12.8 \times 3.7\text{mm}$ $15\mu\text{m}$ filter inside 20 drops $\approx 1\text{ml}$
 Roller clamp
 Latex flash ball
 21G \times 1 1/2" needle connect to luer silp(double wings)
 150cm tubing in total length
 EO sterilization
 Available material :PVC,DEHP FREE PVC
 Available package:Paper and plastic pouch ,PE pouch
 Reference standaed:ISO 8536-4:2010

Item code:JF002
 Specification :
 Air vented spike with cap
 48mm drip with $\text{Ø}12.8 \times 3.7\text{mm}$ $15\mu\text{m}$ filter inside 20 drops $\approx 1\text{ml}$
 Roller clamp
 Latex flash ball
 21G \times 1 1/2" needle connect to luer silp
 150cm tubing in total length
 EO sterilization
 Available material :PVC,DEHP FREE PVC
 Available package:Paper and plastic pouch ,PE pouch
 Reference standaed:ISO 8536-4:2010



Item code:JF003
 Specification :
 Air vented spike with cap
 48mm drip with $\text{Ø}12.8 \times 3.7\text{mm}$ $15\mu\text{m}$ filter inside 20 drops $\approx 1\text{ml}$
 Roller clamp
 Latex flash ball
 21G \times 1 1/2" needle connect to luer lock
 150cm tubing in total length
 EO sterilization
 Available material :PVC,DEHP FREE PVC
 Available package:Paper and plastic pouch ,PE pouch
 Reference standaed:ISO 8536-4:2010



Item code:JF004
 Specification :
 Spike with cap without air-vent
 48mm drip with $\text{Ø}12.8 \times 3.7\text{mm}$ $15\mu\text{m}$ filter inside 20 drops $\approx 1\text{ml}$
 Roller clamp
 Latex flash ball
 21G \times 1 1/2" needle connect to luer lock
 150cm tubing in total length
 EO sterilization
 Available material :PVC,DEHP FREE PVC
 Available package:Paper and plastic pouch ,PE pouch
 Reference standaed:ISO 8536-4:2010



Item code:JF007
 Specification :
 Spike with cap without air-vented
 48mm drip with $\text{Ø}12.8 \times 3.7\text{mm}$ $15\mu\text{m}$ filter inside 20 drops $\approx 1\text{ml}$
 Roller clamp
 Latex flash ball
 21G \times 1 1/2" needle connect to luer slip
 150cm tubing in total length
 EO sterilization
 Available material :PVC,DEHP FREE PVC
 Available package:Paper and plastic pouch ,PE pouch
 Reference standaed:ISO 8536-4:2010

Item code:JF005
 Specification :
 Spike with cap without air-vent
 48mm drip with $\text{Ø}12.8 \times 3.7\text{mm}$ 20 drops $\approx 1\text{ml}$ without filter
 Roller clamp
 Latex flash ball
 21G \times 1 1/2" needle connect to luer lock
 150cm tubing in total length
 EO sterilization
 Available material :PVC,DEHP FREE PVC
 Available package:Paper and plastic pouch ,PE pouch
 Reference standaed:ISO 8536-4:2010



Item code:JF008
 Specification :
 Spike with cap without air-vent
 48mm drip with $\text{Ø}12.8 \times 3.7\text{mm}$ 20 drops $\approx 1\text{ml}$ without filter
 Roller clamp
 Latex flash ball
 21G \times 1 1/2" needle connect to luer slip
 150cm tubing in total length
 EO sterilization
 Available material :PVC,DEHP FREE PVC
 Available package:Paper and plastic pouch ,PE pouch
 Reference standaed:ISO 8536-4:2010



Item code:JF006
 Specification :
 Air vented spike with cap
 48mm drip with $\text{Ø}12.8 \times 3.7\text{mm}$ 20 drops $\approx 1\text{ml}$ without filter
 Roller clamp
 Latex flash ball
 21G \times 1 1/2" needle connect to luer slip
 150cm tubing in total length
 EO sterilization
 Available material :PVC,DEHP FREE PVC
 Available package:Paper and plastic pouch ,PE pouch
 Reference standaed:ISO 8536-4:2010



Item code:JF009
 Specification :
 Air vented spike with cap
 48mm drip with $\text{Ø}12.8 \times 3.7\text{mm}$ $15\mu\text{m}$ filter inside 20 drops $\approx 1\text{ml}$
 Roller clamp
 21G \times 1 1/2" needle connect to luer lock
 150cm tubing in total length
 EO sterilization
 Available material :PVC,DEHP FREE PVC
 Available package:Paper and plastic pouch ,PE pouch
 Reference standaed:ISO 8536-4:2010



Item code:JF010
 Specification :
 Air vented spike with cap
 48mm drip with Ø12.8*3.7mm 15µm filter inside 20 drops ≈ 1ml
 Roller clamp
 Y-connector
 Luer lock with cap
 150cm tubing in total length
 EO sterilization
 Available material :PVC,DEHP FREE PVC
 Available package:Paper and plastic pouch ,PE pouch
 Reference standaed:ISO 8536-4:2010




Item code:JF013
 Specification :
 Air vented spike(double wings) with cap
 48mm drip with Ø12.8*3.7mm 15µm filter inside 20 drops ≈ 1ml
 Roller clamp
 Latex flash ball
 21G×1 1/2" needle connect to luer slip
 150cm tubing in total length
 EO sterilization
 Available material :PVC,DEHP FREE PVC
 Available package:Paper and plastic pouch ,PE pouch
 Reference standaed:ISO 8536-4:2010

Item code:JF011
 Specification :
 Air vented spike with cap
 48mm drip with Ø12.8*3.7mm 15µm filter inside 20 drops ≈ 1ml
 Roller clamp
 Latex flash ball
 Y-connector
 21G×1 1/2" needle connect to luer silp
 150cm tubing in total length
 EO sterilization
 Available material :PVC,DEHP FREE PVC
 Available package:Paper and plastic pouch ,PE pouch
 Reference standaed:ISO 8536-4:2010



Item code:JF014
 Specification :
 Air vented spike(double wings) with cap
 48mm drip with Ø12.8*3.7mm 15µm filter inside 20 drops ≈ 1ml
 Roller clamp
 Latex flash ball
 21G×1 1/2" needle connect to luer lock
 150cm tubing in total length
 EO sterilization
 Available material :PVC,DEHP FREE PVC
 Available package:Paper and plastic pouch ,PE pouch
 Reference standaed:ISO 8536-4:2010



Item code:JF012
 Specification :
 Air vented spike with cap
 48mm drip with Ø12.8*3.7mm 15µm filter inside 20 drops ≈ 1ml
 Roller clamp with wing
 Latex flash ball
 21G×1 1/2" needle connect to luer lock
 150cm tubing in total length
 EO sterilization
 Available material :PVC,DEHP FREE PVC
 Available package:Paper and plastic pouch ,PE pouch
 Reference standaed:ISO 8536-4:2010



Item code:JF015
 Specification :
 Air vented spike(double wings) with cap
 48mm drip with Ø12.8*3.7mm 15µm filter inside 20 drops ≈ 1ml
 Roller clamp
 Latex flash ball
 21G×1 1/2" needle connect to luer slipwith wings
 150cm tubing in total length
 EO sterilization
 Available material :PVC,DEHP FREE PVC
 Available package:Paper and plastic pouch ,PE pouch
 Reference standaed:ISO 8536-4:2010



Item code:JF016
 Specification :
 Air vented spike with cap
 48mm drip with $\text{Ø}12.8 \times 3.7\text{mm}$ $15\mu\text{m}$ filter inside 20 drops $\approx 1\text{ml}$
 Roller clamp with wing
 Latex flash ball
 21G \times 1 1/2" needle connect to luer slip
 150cm tubing in total length
 EO sterilization
 Available material :PVC,DEHP FREE PVC
 Available package:Paper and plastic pouch ,PE pouch
 Reference standaed:ISO 8536-4:2010

Item code:JF017
 Specification :
 Air vented spike with cap
 48mm drip with $\text{Ø}12.8 \times 3.7\text{mm}$ $15\mu\text{m}$ filter inside 20 drops $\approx 1\text{ml}$
 Roller clamp
 Latex flash ball
 22G scalp vein set connect to luer slip
 150cm tubing in total length
 EO sterilization
 Available material :PVC,DEHP FREE PVC
 Available package:Paper and plastic pouch ,PE pouch
 Reference standaed:ISO 8536-4:2010



Item code:JF018
 Specification :
 Air vented spike with cap
 48mm drip with $\text{Ø}12.8 \times 3.7\text{mm}$ $15\mu\text{m}$ filter inside 20 drops $\approx 1\text{ml}$
 Roller clamp
 Latex flash ball
 21G scalp vein set connect to luer slip
 150cm tubing in total length
 EO sterilization
 Available material :PVC,DEHP FREE PVC
 Available package:Paper and plastic pouch ,PE pouch
 Reference standaed:ISO 8536-4:2010



Item code:JF019
 Specification :
 spike with cap without air-vent
 48mm drip with $\text{Ø}12.8 \times 3.7\text{mm}$ 20 drops $\approx 1\text{ml}$ without filter
 Roller clamp
 Latex
 22G scalp vein set connect to luer slip
 150cm tubing in total length
 EO sterilization
 Available material :PVC,DEHP FREE PVC
 Available package:Paper and plastic pouch ,PE pouch
 Reference standaed:ISO 8536-4:2010

Item code:JF020
 Specification :
 Spike with cap without air-vent
 48mm drip with $\text{Ø}12.8 \times 3.7\text{mm}$ 20 drops $\approx 1\text{ml}$ without filter
 Roller clamp
 Latex
 23G scalp vein set connect to luer slip
 150cm tubing in total length
 EO sterilization
 Available material :PVC,DEHP FREE PVC
 Available package:Paper and plastic pouch ,PE pouch
 Reference standaed:ISO 8536-4:2010



Item code:JF019
 Specification :
 spike with cap without air-vent
 48mm drip with $\text{Ø}12.8 \times 3.7\text{mm}$ 20 drops $\approx 1\text{ml}$ without filter
 Roller clamp
 Latex
 21G scalp vein set connect to luer slip with wings
 150cm tubing in total length
 EO sterilization
 Available material :PVC,DEHP FREE PVC
 Available package:Paper and plastic pouch ,PE pouch
 Reference standaed:ISO 8536-4:2010

Blood Transfusion Set

Item code:JFS01
 Specification :
 Air vented spike with cap
 95mm drip with Ø17.3*65mm 200µm filter inside 20 drops ≈ 1ml
 Roller clamp
 Latex flash ball
 21G×1 1/2" needle connect to luer silp(double wings)
 150cm tubing in total length
 EO sterilization
 Available material :PVC,DEHP FREE PVC
 Available package:Paper and plastic pouch ,PE pouch
 Reference standaed:ISO 8536-4:2010



Item code:JFS02
 Specification :
 Air vented spike with cap
 95mm drip with Ø17.3*65mm 200µm filter inside 20 drops ≈ 1ml
 Roller clamp
 Latex flash ball
 19G needle connect to luer lock
 150cm tubing in total length
 EO sterilization
 Available material :PVC,DEHP FREE PVC
 Available package:Paper and plastic pouch ,PE pouch
 Reference standaed:ISO 8536-4:2010

Item code:JFS03
 Specification :
 Air vented spike with cap
 95mm drip with Ø17.3*65mm 200µm filter inside 20 drops ≈ 1ml
 Roller clamp
 19G needle connect to luer lock
 150cm tubing in total length
 EO sterilization
 Available material :PVC,DEHP FREE PVC
 Available package:Paper and plastic pouch ,PE pouch
 Reference standaed:ISO 8536-4:2010



Stopcock

Handle Color:
 Blue,red and white are available
 Two female luer locks and one male luer lock
 Left and right hands are available
 No leakage under pressure:100PSI
 Available package:
 Blister package ,paper and plastic pouch
 Reference standard:ISO 856-10:2004



Extension Line

Specification :
 1.Assemble stopcock with tube
 2.Tube dimenson:3*4.0mm,could be adjusted as request
 3.Tube length could be adjusted ad request
 Avaliable package :paper and plastic pouch

Heparin Cap

Specification :
 Universal 6% taper ,compatible with standard product
 mainly used with I.V.catheter together



Scalp Vein Set

Item code: JF V001

Specification :

Needle point is exquisite, sharp, easy puncture, less pain and little damage

Cannula with thin wall and large inner dimension, smooth flow

Colorful wings for needle size identification

Different sizes meet different requirements

PE /paper individual package

Needle sizes adjust based on customers' requests



IV cannula /IV catheter

IV cannula /IV catheter pen-like type

1.Type :pen-like shape type

2.Needle size :14G -26G


3.Packing :Cap ,paper or blister package



IV cannula /IV catheter with wings
 1.Type :Butterfly wings type
 2.Needle size:14G -26G
 3.Packing :hard blister package



IV cannula /IV catheter with injection port and with wings
 1.Type :with injection port and with wings
 2.Needle size:14G -26G
 3.Packing :hard blister package



IV catheter with small wings

- 1:Type with small wings
- 2:Needle size:14G -26G
- 3:Packing :paper or blister package



Blunt needle

- 1.Size: 21G ,22G ,23G ,25G ,27G
- 2.Packing :paper package



Spinal needle

- 1:Spinal needle used in subarachnoid puncture for spinal anesthesia
- 2:Spinal needle transport hub
- 3:Easy to observe the cerebrospinal fluid discharge,needle core hub with lock can keep the surface of needle bevel joint well
- 4:The needle core whti pin allows a good joint well
- 5:Size :18-26G



Burette infusion set

