

Medical Tape Series



NINGBO HAOSHU HAORUN MEDICAL DRESSING CO.,LTD

ADD: NO1340 LIYUAN ROAD(N) HANSHU NINGBO ZHEJIANG CHINA 315000

Tel:0086-574-87706193

Fax:0086-574-87706193

Mobile:0086-13566542343

Email:sales@haorunmed.com

Website:www.haorunmed.com



NINGBO HAOSHU HAORUN MEDICAL DRESSING CO.,LTD

Bandages for supermarket



Fabric

Size: 10x5cm

Package: 10pcsx100boxes



Waterproof

Size: 45x9mm

Package: 10pcsx100boxes



Fabric

Size: 76x38mm&50x45mm

Package: 20pcsx100boxes



Waterproof

Size: 72x19mm

Package: 30pcsx100boxes

Bandages for supermarket



Waterproof

Size: 76x51mm

Package: 10pcsx100boxes



Waterproof

Size: Dia22mm

Package: 50pcsx100boxes



Waterproof

Size: 76x25mm

Package: 20pcsx100boxes



Waterproof

Size: 56x19mm&76x19mm&76x25mm

Package: 40pcsx100boxes



First Aid Kit

2 Antiseptic cleaning wipes

4 non-woven pad

2 Fabric Bandages 76x19mm

4 Fabric Bandages 56x19mm

Bandages



Material: PE/PVC/EVA/PU Transparent
 Package: 20/50/100pcs/box 300/100/100boxes/carton



Material: Cartoon PE/PVC
 Package: 100/200/300boxes/carton



Material: Fabric/Elastic Fabric/Cotton
 Package: 20/50/100pcs/box 300/100/100boxes/carton

Bandages



Size: 4x10cm/6x10cm/8x10cm customizable
 Material: Plastic
 Package: 10pcs/box 100boxes/carton



Size: 4x10cm/6x10cm/8x10cm customizable
 Material: Fabric
 Package: 10pcs/box 100boxes/carton



Size: 4x10cm/6x10cm/8x10cm customizable
 Material: Non-woven
 Package: 10pcs/box 100boxes/carton



Size: Spot22mm/40X10mm/72X19mm/76X19mm customizable
 Material: PE/PVC/EVA/PU
 Package: 50/75/100pcs/box 50/100/200boxes/carton



Size: Spot22mm/40X10mm/72X19mm/76X19mm customizable
 Material: Fabric/Elastic Fabric/Cotton
 Package: 50/75/100pcs/box 50/100/200boxes/carton

Bandage



Bandage



STERILE WOUND DRESSING



Item Name	Size	MOQ	Package	carton size	Item Name	Size	MOQ	Package	carton size
Wound Dressing	6cmx7cm	30000pcs/size	50pcs*24boxes/carton	30*26.5*28cm	Wound Dressing	9cmx25cm	30000pcs/size	50pcs*18boxes/carton	62.5*55.5*32cm
	6cmx8cm		50pcs*24boxes/carton	30*26.5*28cm		9cmx30cm		50pcs*12boxes/carton	37.5*37*29cm
	6cmx10cm		50pcs*24boxes/carton	39.5*30.5*29cm		9cmx35cm		50pcs*12boxes/carton	42*37.5*29cm
	8cmx10cm		50pcs*24boxes/carton	55.5*32.5*32cm		10cmx10cm		50pcs*24boxes/carton	55.5*32.5*32cm
	8cmx15cm		50pcs*24boxes/carton	55.5*42.5*32cm		10cmx12cm		50pcs*24boxes/carton	55.5*42.5*32cm
	8cmx20cm		50pcs*24boxes/carton	55.5*52.5*32cm		10cmx15cm		50pcs*24boxes/carton	55.5*42.5*32cm
	9cmx10cm		50pcs*24boxes/carton	55.5*32.5*32cm		10cmx20cm		50pcs*24boxes/carton	55.5*52.5*32cm
	9cmx15cm		50pcs*24boxes/carton	55.5*42.5*32cm		10cmx25cm		50pcs*18boxes/carton	62.5*42*32cm
	9cmx20cm		50pcs*24boxes/carton	55.5*52.5*32cm		10cmx30cm		50pcs*12boxes/carton	55.5*37*32cm
	9cmx35cm		50pcs*12boxes/carton	55.5*42*32cm					

Material: Non-woven/pu

SRERILE WOUND DRESSING



PU/Non woven fixing Tape Roll

Normal Size	Package	Carton size
5cm x 10m	1roll/box,160boxes/ctn	61.5x39x23cm
10cm x 10m	1roll/box,80boxes/ctn	39x31.5x43cm
15cm x 10m	1roll/box,60boxes/ctn	39x31.5x47.5cm
20cm x 10m	1roll/box,40boxes/ctn	39x31.5x42.5cm

We can do as customer request

MEDICAL ADHESIVE TAPE



Non-Woven Tape(Skin/White Non-Woven)			
Size	Package	Carton size	Moq
1.25cm x 5m	24rolls/box,30boxes/ctn	25x24.5x28cm	1500Boxes/Each Size
2.5cm x 5m	12rolls/box,30boxes/ctn		
5cm x 5m	6rolls/box,30boxes/ctn		
7.5cm x 5m	6rolls/box,30boxes/ctn		
10cm x 5m	6rolls/box,30boxes/ctn		



Transparent Pe Tape(Transparent)			
Size	Package	Carton size	Moq
1.25cm x 5m	24rolls/box,30boxes/ctn	27x26.5x28cm	1500Boxes/Each Size
2.5cm x 5m	12rolls/box,30boxes/ctn		
5cm x 5m	6rolls/box,30boxes/ctn		
7.5cm x 5m	6rolls/box,30boxes/ctn		
10cm x 5m	6rolls/box,30boxes/ctn		



Silk Tape Sawteeth(Silk)			
Size	Package	Carton size	Moq
1.25cm x 5m	24rolls/box,30boxes/ctn	27x26.5x28cm	1500Boxes/Each Size
2.5cm x 5m	12rolls/box,30boxes/ctn		
5cm x 5m	6rolls/box,30boxes/ctn		
7.5cm x 5m	6rolls/box,30boxes/ctn		
10cm x 5m	6rolls/box,30boxes/ctn		



Zinc Oxide Adhesive Tape(Skin/White Cotton)			
Size	Package	Carton size	Moq
1.25cm x 5m	24rolls/box,30boxes/ctn	32x32x30cm	1500Boxes/Each Size
2.5cm x 5m	12rolls/box,30boxes/ctn		
5cm x 5m	6rolls/box,30boxes/ctn		
7.5cm x 5m	6rolls/box,30boxes/ctn		
10cm x 5m	6rolls/box,30boxes/ctn		

MEDICAL ADHESIVE TAPE



Zinc Oxide Adhesive Tape(Skin/White Cotton)			
Size	Package	Carton size	Moq
1"	1/2" x 4 rolls	12" x 10yds	1500Boxes/Each Size
	1" x 3 rolls		
	2" x 2 rolls		



Zinc Oxide Perforated Plaster			
Size	Package	Carton size	Moq
10cm x 5m	1roll/box,80boxes/ctn	45.5x29x25cm	3000Rolls/Each Size
10cm x 5m	10rolls/box,80boxes/ctn	36.5x36.5x38.5cm	
18cm x 5m	1roll/box,80boxes/ctn	51.5x41.5x21.5cm	
18cm x 5m	10rolls/box,10boxes/ctn	53x51.5x22cm	



Zinc Oxide Tinplate Plaster(Skin/White Cotton)			
Size	Package	Carton size	Moq
1.25cm x 5m	12rolls/box,60boxes/ctn	54x48x27cm	1500Boxes/Each Size
2.5cm x 5m	12Rolls/box,30boxes/ctn	36.5x36.5x38.5cm	
5cm x 5m	6Rolls/box,30boxes/ctn	36.5x36.5x33.5cm	
7.5cm x 5m	6Rolls/box,20boxes/ctn	36.5x24.5x45.5cm	



Zinc Oxide Plaster with plastic cover(Skin/White Cotton)			
Size	Package	Carton size	Moq
1.25cm x 5m	24rolls/box,30boxes/ctn	36.5x36x42.5cm	1500Boxes/Each Size
2.5cm x 5m	12rolls/box,30boxes/ctn	36.5x36x34.5cm	
5cm x 5m	6rolls/box,30boxes/ctn	36.5x36x31cm	
7.5cm x 5m	6rolls/box,30boxes/ctn	36.5x36x44.5cm	
10cm x 5m	6rolls/box,30boxes/ctn	59.5x36x34.5cm	

MEDICAL SPORT TAPE



EAB Tape(Cotton)		
Size	Package	Carton Size
5cm x 4.5m	144rolls	46x31x36cm
7.5cm x 4.5m	96rolls	46x31x36cm
10cm x 4.5m	96rolls	46x31x46cm
15cm x 4.5m	48rolls	46x31x36cm
5cm x 3m	144rolls	38.5x25.5x36cm
7.5cm x 3m	96rolls	38.5x25.5x36cm
10cm x 3m	96rolls	38.5x25.5x46cm
15cm x 3m	48rolls	38.5x25.5x36cm

Sport Tape(Cotton)		
Size	Package	Moq
2.5cm x 5m	12Rolls/box,30boxes/ctn	1500Boxes/Each Size
5cm x 5m	6Rolls/box,30boxes/ctn	
3.8cm x 5m		



Non-Woven Cohesive Bandage(printing)		
Size	Package	Moq
2.5cm x 4.5m	24Rolls/box,24boxes/ctn	1500Boxes/Each Size
5cm x 4.5m	12Rolls/box,524boxes/ctn	
7.5cm x 4.5m	12Roll/box,16boxes/ctn	
10cm x 4.5m	12Roll/box,8boxes/ctn	
15cm x 4.5m		

MEDICAL SPORT TAPE



B. Strong stickines with acrylic glue (latex free), more than 3 days stay in place and no residue & irritation to skin

A. High elasticity of 170%~180%




D. Waterproof and no release with sweating, can bath after soprts

C. High quality breathable cotton with spandex, similar to the thickness and elasticity of human skin, provide the best comfort



Kinesiology Tape(95% Cotton+5%Spandex)			
Size	Package	Carton size	MOQ
2.5cm x 5m	12rolls/box,20boxes/ctn	39x30.5x35cm	5Cartons/Each Size
3.8cm x 5m	8rolls/box,20boxes/ctn	39x30.5x35cm	
5cm x 5m	6rolls/box,20boxes/ctn	39x30.5x35cm	
7.5cm x 5m	4rolls/box,20boxes/ctn	39x30.5x35cm	
10cm x 5m	12rolls/box,6boxes/ctn	48.5x33.5x36cm	
We can do as customer request.			

COOLING PATCH



COOLING GEL PACTH

Size: 12x5cm/11x4cm
 Material: Non-woven+Gel+Film
 Package: 1patchx3bagsx240boxes/ctn
 MOQ: 90000bags



ANTISEPTIC CLEAN WIPE



Alcohol Pad(Non-Woven Fabric+Alcohol)		
Size	Package	Moq
6cm x 3cm	100pcsx100boxes/carton	500000PCS
	200pcsx50boxes/carton	

Povidone-Iodine Prep Pad(Non-Woven Fabric+Povidone-Iodine prep)		
Size	Package	Moq
6cm x 3cm	100pcsx100boxes/carton	500000PCS



Mouth Wash/Hand Sanitizer




Size:100ml/250ml
 500ml



Size:60ml

Capsicum Plaster



Capsicum Plaster

6x10cm	50bags/box	12boxes/CTN	CTN Size:32x21.5x17cm
10x18cm	50bags/box	12boxes/CTN	CTN Size:32x30x25cm
12x18cm	24bags/box	50boxes/CTN	CTN Size:57x36x25cm

Cotton swab



cotton swab

Material:plasticod paper wood
Size:200pcs x12box x10shrink

cotton swab

Material:plastic/paper/wood
Size:300pcs x12box x30shrink

cotton swab

Material:plastic/paper/wood
Size:100pcs x per bag, 300bags x per carton

cotton swab

Material:plastic/paper/wood
Size:100pcs x240boxes